



WellCare™ ATH[®] Systems

Automatic Tightening Band Ligation for Hemorrhoids



**Automatic
Quantitative
Self-tightening**

The **most effective** clinical solution for grade I and II and select patients with grade III internal hemorrhoidal disease¹

Patented² WellCare™ ATH[®] System is an **automatic, convenient, and efficient** non-surgical treatment method for all indicated hemorrhoids in one procedure

Reference


¹ The American Society of Colon and Rectal Surgeons Clinical Practice Guidelines for the Management of Hemorrhoids , 2018

² International Patent No. PCT/CN2017089297

WellCare™ ATH® Systems

Automatic Tightening Band Ligation for Hemorrhoids

Patented WellCare™ ATH® System is an automatic, convenient, and efficient non-surgical treatment method for grade I and II and select patients with grade III internal hemorrhoidal disease. Its design allows treatment of all indicated hemorrhoids in one procedure.



Ligator Barrel

Replaceable ligator barrel with transparent and visible suction chamber improve the success of ligation procedure for hemorrhoids and meet various clinical requirements.

Pusher

Unique pusher designs allow elastic thread to maintain sustainable and quantitative force which ensure standard of care.

Self-tightening Trigger

Patented designs allows tightening and ligation process in one procedure by one hand operation.

Order Information



Model	Description	Unit of Measure	Minimum Oder Quantity(Set)
ATET-02B-3H	1 Ligator Body+3 Barrels/Auto Tightening and Manual cutting	Set	100
ATET-02B-4H	1 Ligator Body+4 Barrels/Auto Tightening and Manual cutting	Set	100
ATET-02C-3H	1 Ligator Body+3 Barrels/Manual Tightening and Manual cutting	Set	100
ATET-02C-4H	1 Ligator Body+4 Barrels/Manual Tightening and Manual cutting	Set	100
ATET-03C-3H	1 Ligator +3 Elastic String Reloaders/Manual Tightening and Manual cutting	Set	100
ATET-03C-4H	1 Ligator +4 Elastic String Reloaders/Manual Tightening and Manual cutting	Set	100
ATET-02D-4H	1 Ligator Body+4 Barrels/Auto Tightening and Auto Cutting	Set	100
1H	Single ligator Barrel	Pcs	200

www.wellcaremedtech.com

Well Care Medical Technology Co., Ltd.

ADD: No. 11# The No.1 Road Fozuling Of East Lake High-Tech Development Zone, Wuhan, Hubei, China
China Wuhan Optical Valley

TEL/FAX: +86-027-65522456

EMAIL: info@wellcaremedtech.com