



Products Catalogue

A Medical Reference Guide



WellCare™ ANO System

Advanced Anoscopy Solutions


- Wellcare™ANO offers a wide range of disposable self-illuminating anoscopes for use in the diagnosis and treatment of conditions of the anal canal and lower rectum.
-
- Wellcare™ANO System consists of Seven main pipeline : Wellcare™ANO Light, Wellcare™ANO,Wellcare™ANO Truncated WellCare™ ANO Truncated (AS-1) , WellCare™ ANO Integrate,WellCare™ ANO, WellCare™ ANO Beveled.

WellCare™ ANO Light

Led Anoscope

Wellcare™ANO Light is a disposable, self-illuminating anoscope that functions as an all-in-one integrated system for diagnosis, treatment, and follow-up. It enhances examination efficiency, reduces the need for multiple instruments, supports modern proctology clinical workflows, and improves procedural accuracy and convenience.

Removable Self –illuminated LED light
Provides the highest lighting quality to ensure the success of operation.





Scaled LED Integrated comfort sleeve
Provides the clinician with the flexibility to posit the Proctoscopy at any angel. It makes a convenient clinical application.

Ergonomic Trocar design
Allows for enhancing patient comfort.

WellCare™ ANO

WellCare™ ANO Truncated (AS-1) ►

WellCare™ ANO Truncated features a large-diameter, short length design that ensures smooth insertion while significantly expanding the surgical field, enhancing the convenience of examination and surgical procedures.




WellCare™ ANO Beveled

WellCare™ ANO Beveled includes a speculum and a sleeve, designed to guide instruments and protect surrounding tissue during examinations.

◀ WellCare™ ANO Integrate

WellCare™ ANO Integrate consists of an trocar and a speculum, and is additionally equipped with a sleeve for guided instrument insertion, helping to protect surrounding tissue during procedures. In proctological diagnosis and treatment, it enables precise instrument positioning and improves procedural accuracy.



◀ WellCare™ ANO

WellCare™ ANO is a basic disposable anoscope composed of a speculum and an trocar. Its simple and reliable structure allows for easy operation, making it suitable for routine anorectal examinations and basic diagnostic procedures.


WellCare™ ANO

Led Anoscope for Laser treatment

Specially designed for laser treatment, Wellcare™ ANO features a stable structure that ensures precision during operations. It offers optimized access and excellent visibility throughout therapeutic procedures, while the integrated self-illuminating design enhances visualization without the need for external light sources.

Removable Self –illuminated LED light

Provides the highest lighting quality to ensure the success of operation.



Lateral Slot and Scaled design

Designed for laser treatment, it offers unrestricted angular flexibility to position the anoscope, ensuring optimal precision and ease during treatment.

Ergonomic Trocar design

Allows for enhancing patient comfort.


WellCare™ ANO Truncated

Led Anoscope

Wellcare™ANO Truncated is a disposable anoscope designed for basic anorectal examination and routine anoscopy procedures. Its standard design, compatible with LED light sources, ensures clear visualization—making it a reliable choice for primary and general clinical settings.

Removable Self –illuminated LED light

Provides the highest lighting quality to ensure the success of operation.



Lateral Slot and Scaled design

Standard design for basic diagnostic examination and routine anoscopy procedures.

Ergonomic Trocar design

Allows for enhancing patient comfort.



www.wellomedical.com

Well Care Medical Technology Co.,Ltd.


ADD:No.11#The No.1 Road Fozuling Of
East Lake High-Tech Development Zone,
Wuhan,Hubei,China
China Wuhan Opical Valley

TEL/FAX:+86-027-65522456

EMAIL:info@wellomedical.com

EXCELLENCE IN PELVIC DISEASE THERAPY

Order Information

PRODUCT PICTURE	MODEL	DESCRIPTION	UNIT OF MEASURE	QUANTITY
	WellCare-ANO-Light (AS-20)	An all-in-one integrated anoscope system for diagnosis, treatment, and follow-up.	Box	100
	WellCare-ANO-Truncated (AS-30)	Suitable for routine anorectal examination.	Box	100
	WellCare-ANO (AS-40)	Specially designed for laser treatment.	Box	100
	WellCare-ANO-Truncated (AS-1)	Large-diameter and short-length design facilitates surgical procedures.	Box	100
	WellCare-ANO-Integrate (AS-2)	Equipped with a guide sleeve for safe instrument insertion and tissue protection.	Box	100
	WellCare-ANO-Beveled (AS-3)	Simplified design with sleeve for instrument guidance and tissue protection.	Box	100
	WellCare-ANO	Simple, reliable structure for operation in routine anorectal examinations.	Box	100